


# TruAir™ Scope Storage Wall Cabinet

<b>Model Numbers:</b>	M1201HG-HEPA M1201HG-HEPA-RIGID-HYB	<b>Ext. Dimensions:</b>	48-3/8"H X 19-1/4"W X 11-3/4"D
<b>Color:</b>	Bright White	<b>Construction:</b>	HDPE Plastic Locking Hinged Glass Doors
<b>Cabinet Filter Efficiency:</b>	Removes 99.97% of Particles (up to 0.3µm)	<b>Air Supply:</b>	HEPA Filter Fan (Cabinet)
<b>Capacity:</b>	4-8 Flexible or Rigid Small Diameter Endoscopes	<b>Electrical Supply:</b>	120 V Wall Plug-3 A (6' Cord Exiting Top Back of Cabinet)
		<b>Options:</b>	TruHEPA Air System Slope Top


## Measurements

Item 1  
PLASTIC CABINET  
M1201HG-HEPA Capsa HALF SIZE WALL SCOPE CABINET  
with Hinged See-through Door



# TruAir™ Scope Storage Wall Cabinet

## Measurements



Item  
**PLASTIC CABINET**  
**M1201HG-HEPA-RIGID-HYB Capsa HALF SIZE WALL SCOPE CABINET**  
 with Hinged See-through Door