



# MEDICATION MANAGEMENT AND ORGANIZATION

- Medication Carts
- Medical & Procedure Carts
- Automated Dispensing
- Computing & Telehealth/Virtual Care
- Inventory Management

1 This is Capsa Healthcare

A healthcare market leader for over 60 years through inventive product design and commitment to customer support.

### Table of Contents

Avalo Series Medication Carts	5
Avalo Woodblend Medication Carts	9
M-Series Medication Carts	13
Vintage Encore Medication Carts	19
Specialty Carts	23
NexsysADC™ Automated Dispensing Technology	25
Avalo Series Emergency & Medical Carts	29
Computing & Telehealth/Virtual Care	31
Kirby Lester Inventory Management	35
Warranty, Service, & Support	36





Organized storage. Improved productivity. Innovative technology. Capsa Healthcare provides the right medication management products to help you deliver better care in assisted living, skilled nursing, and long-term care settings.

The most comprehensive product suite in the industry is all right here - Medication Carts and Cabinets, Computing Workstations, Medical Carts, Advanced Medication Management Systems, and Inventory Solutions. Capsa precisely delivers what your facility requires!



*Over 99% of Capsa products are delivered on-time to our customers' sites.*

*Over 99% of Capsa products pass our rigorous 55-point QC process the first time.*

### Capsa Healthcare Overview

- Headquartered in Columbus, Ohio, with additional sales, management, and manufacturing operations in Portland, Oregon, and Waalwijk, The Netherlands
- Product installations in over 70 countries and thousands of healthcare facilities worldwide
- Our comprehensive product suite ensures medication administration accuracy and control from the pharmacy to the patient, and greatly improves the efficiency of the nursing staff







# AVALO<sup>®</sup> SERIES MEDICATION CARTS





# Avalo® Series

## Medication Carts

**Avalo Series** medication carts excel in areas that make the difference for nurses, like easy usage and access, a wide range of drawer configurations, enormous capacity, and smooth drawer action. For the pharmacy, Avalo yields consistently high user ratings, proven lock systems, and the ability to refurbish panels for lasting value.

### The Avalo Difference

- Modular design for simple component replacement and refurbishing
- Support virtually all medication delivery systems including punch cards, auto-packaging, unit dose box
- High-density panels prevent denting or corrosion
- 3", 6", and 10" seamless drawers
- Unique integrated drawer divider system
- Full selection of accessories and lock systems
- Maximize work space with several computing mount options for EMR/EHR

◀ *Avalo Medication Carts are long-lasting and a superior value.*



**PCL punch card cart**  
*shown with optional accessories*



**UDS automated packaging cart**  
*shown with optional accessories*

# Punch Card Models

## Model PCS

43.75"(h) 24"(d) 31.5"(w)



## Model PCL

43.75"(h) 24"(d) 38"(w)



## Model PCXL

43.75"(h) 24"(d) 48.5"(w)



Height is floor to work surface. Standard drawers are 3"(d) and 10"(d). Accommodates full and half card sizes, plus liquids & supplies.



Offers a practical capacity of 350 punch cards plus room for liquids and supplies. Maximum capacity 450 punch cards.



Offers a practical capacity of 450 punch cards plus room for liquids and supplies. Maximum capacity 600 punch cards.



Offers a practical capacity of 600 punch cards plus room for liquids and supplies. Maximum capacity 750 punch cards.



Avalo accommodates virtually all popular compliance packaging systems, such as Parata®/TCGRx®, Omnicell®, Dispill®, SynMed® and Opus®

## Inquire about our Treatment Carts

- Large capacity & flexible organization
- Expansive work area
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panel
- 10-year limited warranty

See p. 23 for more



# Auto-Packaging Models

## Model UDS

43.75"(h) 24"(d) 31.5"(w)



## Model UDL

43.75"(h) 24"(d) 38"(w)



## Model UDXL

43.75"(h) 24"(d) 48.5"(w)




Height is floor to work surface. Standard drawers are 3"(d) and 10"(d). Accommodates most medication pass methods.

## Auto-Packaging Drawer Configuration Options

UDL drawers shown. Ask about VersaBin lids and dividers, label holders and alternative bin sizes.

### VersaBin™ Complete



### VersaBin™ Transfer



### VersaBin™/Fixed Divider System\*



### Fixed Divider System



### Drawer Flex System



Ask about optional bin dividers, lids, label holders, and alternative bin sizes.

\*VersaBin is a customizable patient medication organization system, with optional transfer trays and totes





# AVALO<sup>®</sup> WOODBLEND MEDICATION CARTS

## Medication Carts

The **Avalo Woodblend Series** medication cart line features enhanced styling and improved modularity to facilitate component replacement and refurbishing. Modern healthcare providers will appreciate the superior organization of medication administration systems and clinical supplies. Avalo's long-term value provides the pharmacy a lower cost of ownership and improved cost management.

### The Avalo Woodblend Difference

- Modular design for simple component replacement and refurbishing
- High-density panels maintain like-new appearance and offer enhanced durability
- Wide array of configurations available featuring 3", 6", and 10" seamless drawers
- Full selection of computer arm and tray systems for EHR/eMAR software, plus multiple accessories and locking options
- Premium casters simplify cart movement by 20%



**PCS punch card cart**

*shown in Mountain Ash with optional accessories*



**UDL automated packaging cart**

*shown in Northern Cherry with optional accessories*

### Easier to Maneuver



### Mobilize EMR/EHR



### Available In 3 Widths



# Security Options

## Lock System Choices

### AutoLock™ Automatic Relocking System



Automatically secures a motionless unlocked medication cart with a user-defined relock time delay from 1-99 minutes. The system also permits the use of up to 5,000 individual user codes.

**See available upgrades below.**

### Keyless Access with Manual Relocking



Simply enter a 4-digit access code, push down the lock lever and the cart is ready for access. To secure the medication cart, push in the lock lever. Clear visual indication of lock status is apparent from a distance.

### Key Access with Manual Relocking (standard)



Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides you with the flexibility to change keys to the medication cart easily. *Standard on Avalo carts.*


## AutoLock™ Upgrades

### Proximity Access



*An advanced proximity card reader can be selected as an additional security access requirement.*

### Audit Tracking



*Authorized users can track cart activity, including user access and relocking for over 2,000 events.*

**AutoLock™** automatically secures a motionless unlocked medication cart with a user-defined relock time.



# Accessories & Options

\* Included on all Avalo Medication Carts.

S & L Models also include Slide-Out Surface. XL Model includes Tracking Caster.

For additional accessories see [CapsaHealthcare.com](http://CapsaHealthcare.com)



Core Removable Lock\*



Waste with Lid\*



Narcotic Storage\*



5" Brake Caster\*



Cart Handle



Slide-Out Surface



Drawer Flex System



Sharps Container



Clear Top Mat



Drawer Dividers



Bulk Storage



Laptop Arm & Tray



Narcotic Drawer (XL only)



Drawer ID Label



Side Organizer




Corner Organizer



Glove Box



Double - Wide Narc



5" Tracking Caster



Wheel Bumpers

## Avalo Series Accent Colors



Taupe



Horizon Blue



Blush Salmon



Extreme Yellow



Lavender



Cantaloupe



Pure Lime

## Woodblend Finishes



Northern Cherry



Mountain Ash



# M-SERIES MEDICATION CARTS



## Medication Carts

Maximum capacity and organization are hallmarks of Capsa Healthcare's M-Series medication carts. Pharmacies appreciate the multiple locking options and proven security, a fast return on investment, and trouble-free usage. Nurses find M-Series to be simple to access medications and supplies, and effortless to guide and open; plus drawers are infinitely customizable. Finally, sleek styling tops off the popular M-Series.

### The M-Series Difference

- Support virtually all medication delivery systems including punch cards, auto-packaging, unit dose box
- Lightweight, superior durability
- 3", 6", and 10" drawers
- Full selection of accessories
- Choice of lock security systems
- Maximize work space with several computing mount options for EMR/EHR



*M-Series Medication Carts offer attractive styling and huge capacity.*

**M3-AP automated packaging cart**

*shown with optional accessories*



**M4-PC punch card cart**

*shown with optional accessories*



## Punch Card Models

### Model M2-PC

43"(h) 26.7"(d) 23.1"(w)



Maximum capacity  
300 punch cards.

### Model M3-PC

43"(h) 26.7"(d) 29.4"(w)



Maximum capacity  
450 punch cards.

### Model M4-PC

43"(h) 26.7"(d) 36.1"(w)



Maximum capacity  
600 punch cards.

### Model M5-PC

43"(h) 26.7"(d) 45.2"(w)



Maximum capacity  
750 punch cards.

### Model M6-PC

43.5"(h) 25.2"(d) 50.5"(w)



Maximum capacity  
900 punch cards.

### Inquire about our Treatment Carts

- Available in M-Series (M3-TC) and Vintage Encore (VE3-TC)
- Large capacity & flexible organization
- Durable drawers, 5 depths
- Full selection of accessories

See p. 23  
for more



Height is floor to work surface. Standard drawers are 3"(d) and 10"(d). Accommodates full and half card sizes, plus liquids & supplies.

## Auto-Packaging Models

### Model M2-AP

43"(h) 26.7"(d) 23.1"(w)



### Model M3-AP

43"(h) 26.7"(d) 29.4"(w)



### Model M4-AP

43"(h) 26.7"(d) 36.1"(w)



### Model M5-AP

43"(h) 26.7"(d) 45.2"(w)




Height is floor to work surface. Drawers are 3"(d) and 10"(d). Accommodates most medication pass methods.


### Auto-Packaging Drawer Configuration Options

M4-AP drawers shown. Ask about Stor-Flex™ bin lids and dividers, label holders and alternative bin sizes.

#### Stor-Flex™ System\*



#### Stor-Flex™ Transfer



#### Stor-Flex™ Bin/Fixed Rail System



#### Fixed Rail System



#### Drawer Flex System



M-Series accommodates virtually all popular compliance packaging systems, such as Parata®/TCGRx®, Omnicell®, Dispill®, SynMed® and Opus®

\***Stor-Flex** is a customizable patient medication organization system, with optional transfer trays and totes


## Stor-Flex™ Auto-Packaging Solutions

Stor-Flex is a simple, logical approach to auto-packaging organization. Storage is in singular bins that can easily accommodate resident discharges, short-term leaves and room changes.

### Stor-Flex™ Features

You decide how resident medications, PRNs, and small items are stored and re-stocked.


- Bin dividers and flexible labeling
- Large resident capacity per cart
- Flexibility to change from punch card to Stor-Flex configuration



Bin modules available in 3 sizes



Large bin perfect for 7-day med cassette systems



Bins load neatly and securely in Stor-Flex Tote System

## M-Series & Encore Series Security Options



### Key Access with Manual Relocking (standard)

Featuring the Best® Core Removable Lock System, this basic yet proven lock system provides you with the flexibility to change keys to the medication cart easily. *Standard on M-Series and Vintage carts.*



### EMS - Entry Management System

The advanced EMS™ system offers unparalleled security in a proven lock system designed specifically for healthcare carts. This technology upgrade makes the M-Series & Vintage Series lock system the most advanced cart security available in long-term care.


- Keyless access with auto-relock
- Permits 250 unique users and supervisors, programmable via keypad or PC
- Optional LockView™ Software provides the platform for secondary security via proximity card and the ability to audit and manage EMS via a PC



## Accessories & Options

\* Included on all M-Series Medication Carts.  
 M3 & M4 Models also include Slide-Out Surface,  
 Waste with Lid. M5 & M6 Models include Tracking Caster,  
 Waste with Lid.

For additional accessories see [CapsaHealthcare.com](http://CapsaHealthcare.com)



Cart Handle\*



Chair Rail Bumpers\*



5" Brake Caster\*



10" Narcotic Storage\*



Core Removable Lock\*



Waste with Lid




Slide-Out Surface



Pill Crusher



Externals Drawer Divider System



Drawer Flex System



Punch Card Exchange Bins



Organizer Tray



Corner Organizer



Sharps Collector



Stor-Flex Divider System



Gloves/Tissue Holder



Laptop Arm & Tray



Narc Drawer (M5 only)



5" Tracking Caster



Double-Wide Narc

### M-Series Accent Colors



Blue-Gray



Platinum



Teal



Cranberry



Reflex Blue



# VINTAGE ENCORE™ MEDICATION CARTS





# Vintage Encore™

## Medication Carts

Vintage Encore is the the senior care industry's benchmark for furniture-like medication carts. A fresh advancement of Capsa's classic Vintage Series, Vintage Encore carts now have enhanced styling and work-saving features. Vintage Encore is ideal for extended care facilities that value modern décor and elegant touches.

### The Vintage Encore Difference

- Support virtually all medication delivery systems including punch cards, auto-packaging, unit dose box
- Maximize work space with several computing mount options for EMR/EHR
- No-bounce extended work surface
- Sleek 1-piece drawer handles
- Soft-close drawers (3", 6", 10")
- Full selection of accessories and locking systems
- Premium casters simplify cart movement by 20%



**VE3-PC punch card cart**

shown in Boston Cherry with optional accessories



**VE4-AP automated packaging cart**

shown in White Ash with optional accessories



# Vintage Encore

## Punch Card Models

**Model VE2-PC**  
44"(h) 25.2"(d) 21.7"(w)



Shown in  
Chatsworth Oak



Maximum capacity  
300 punch cards.

**Model VE3-PC**  
44"(h) 25.2"(d) 28"(w)



Shown in  
Boston Cherry



Maximum capacity  
450 punch cards.

**Model VE4-PC**  
44"(h) 25.2"(d) 34.7"(w)



Shown in  
White Ash



Maximum capacity  
600 punch cards.

**Model VE4x-PC**  
44"(h) 25.2"(d) 41.4"(w)



Shown in  
Chatsworth Oak



Maximum capacity  
750 punch cards.

**Model VE5-PC**  
44"(h) 25.2"(d) 44.7"(w)



Shown in  
White Ash



Maximum capacity  
750 punch cards.

**Model VE6-PC**  
44"(h) 25.2"(d) 50.5"(w)



Shown in  
Antique Walnut



Maximum capacity  
900 punch cards.

Height is floor to work surface. Standard drawers are 3"(d) and 10"(d). Accommodates full and half card sizes, plus liquids & supplies.

### Inquire about our Medication Cabinets

- Wall mount
- Countertop



# Vintage Encore

## Auto-Packaging Models

**Model VE2-AP**  
44"(h) 25.2"(d) 21.7"(w)



Shown in  
Chatsworth Oak

**Model VE3-AP**  
44"(h) 25.2"(d) 28"(w)



Shown in  
Boston Cherry

**Model VE4-AP**  
44"(h) 25.2"(d) 34.7"(w)



Shown in  
Antique Walnut

**Model VE5-AP**  
44"(h) 25.2"(d) 44.7"(w)



Shown in  
White Ash

Height is floor to work surface. Drawers are 3"(d) and 10"(d). Accommodates most medication pass methods.


### Auto-Packaging Drawer Configuration Options

VE4-AP drawers shown. Ask about Stor-Flex™ bin lids and dividers, label holders and alternative bin sizes.

#### Stor-Flex™ System




#### Stor-Flex™ Transfer




#### Stor-Flex™ Bin/Fixed Rail System\*



#### Drawer Flex System



#### Fixed Rail Divider System



\***Stor-Flex** is a customizable patient medication organization system, with optional transfer trays and totes

Vintage Encore accommodates virtually all popular compliance packaging systems, such as Parata®/TCGRx®, Omnicell®, Dispill®, SynMed® and Opus®



# Vintage Encore

## Accessories & Options

\* Included on all Vintage Encore Medication Carts.  
 VE4x-PC, VE5, & VE6 Models also include Tracking Caster.

▶ See p.16 for Security Options



10" Narcotic Box\*



5" Brake Caster\*



Core Removable Lock\*



Sharps Collector



Side-Mount Organizer



Slide-Out Surface



Waste With Lid



Laptop Arm & Tray



Drawer Dividers-see P. 21




Harley Handle



Biohazard Case



Drawer Flex System



5" Tracking Caster



Corner Organizer




Skirted Base (non-mobile)



Narc Drawer (VE5, VE6 only)

### Available Finishes




Boston Cherry



White Ash




Chatsworth Oak



Antique Walnut

### Handle Options



Brushed Nickel



Satin Black

### Cart Top

Durable cart tops that are scratch, gouge and stain resistant





# Specialty Carts

## Treatment Carts

Lightweight performance, simple maneuverability and flexible organization of treatment supplies meet any care facility's storage needs.

- Large capacity & flexible organization of treatment supplies
- Wide selection of accessories
- Expansive work area



**Avalo® Series** (3 Heights)  
 Standard 42.5"(h) 24"(d) 31"(w)  
 Intermediate 39.5"(h) 24"(d) 31"(w)  
 Compact 36"(h) 24"(d) 31"(w)



**Vintage Encore VE3-TC**  
 40.1"(h) 25.2"(d) 28"(w)



**M-Series M3-TC**  
 40.1"(h) 25.2"(d) 28"(w)

## Multi-Dose Medication Carts

Capsa offers three models of medication carts to accommodate most large-format, one-week multi-dose blister packaging systems.

- Available in multiple widths, 21" to 48.5" w
- Deep drawers (12" and 13" h) and unique dividers organize large cards



**M-Series M4-DP**  
 43.5"(h) 25.2"(d) 34.7"(w)



**Vintage Encore VE5-DP**  
 44"(h) 25.2"(d) 44.7"(w)



**Avalo® AVDPL**  
 43.75"(h) 24"(d) 38"(w)



# NEXSYSADC™ AUTOMATED DISPENSING TECHNOLOGY



## Practical Automated Dispensing Technology

NexsysADC™ is the simplest, most secure technology to fully manage the onsite storage and distribution of medications (especially controls), STAT/first doses, and supplies. With NexsysADC, the level of security you need and expect is only a fraction of the price of alternative automated dispensing cabinets. The pharmacy benefits from a faster ROI without compromising safety, security, or operational flexibility. And the facility benefits from having the new Stat or 1st dose on-hand at the exact time it's needed.

### The NexsysADC Difference

- **Medication Security:** Innovative solution to manage the restocking and access of controls and high-value medications.
- **Fast ROI Without Compromise:** Your Capsa Representative will run a detailed analysis to show how NexsysADC can be deployed in more facilities for broader control and higher customer satisfaction, compared to alternative ADCs.
- **Efficient Stat Dose/E-kit Delivery:** Eliminates manual methods like the “tackle box” commonly used in extended care settings. Significantly lowers costs associated with stat deliveries, and improves audit accuracy (no more missing paperwork).
- **Configurable For Each Facility:** Available in a standard cabinet and new countertop unit. Each NexsysADC is custom-built for your specific needs and patient population.
- **Streamlined Replenishment & Inventory:** NexsysADC identifies what medications need to be restocked. The flexible replenishment process saves time, reduces inventory turns, and prevents diversion.

When more capacity is needed, add one or more Auxiliary Cabinets (right) for storing hundreds to thousands of SKUs. Networks to any Main or 4T cabinet.



**NexsysADC** fits your formulary and budget, in main cabinet (left) and 4T countertop. Security and storage for any size care environment:

- Extended care nursing homes
- Critical access hospitals
- Surgery centers
- Group homes
- Rehab and psych environments
- Veterinary clinics
- Dental surgery
- Hospice
- Education simulation






## Guided Access With 'CAM' Technology

NexsysADC relies on guided access with multiple levels of illuminated and on-screen prompts. This ensures the caregiver IDs the exact location of the medication, and picks the right medication for the right patient. Guided access technology also ensures restocking and inventory control are accurately performed.

### Three Secure Ways to Store & Track Controls and High-Value Supplies



**1 Controlled Meds**  
Controlled and high-value medications are safely guarded in CAMs (controlled access modules). Each 3" CAM holds up to 25 smaller SKUs, and deep 6" and 10" CAM Drawers hold up to 12 large-item SKUs.



**2 Routine Meds & Small Supplies**  
Meds and surgery supplies are efficiently managed in dividable, locked open-matrix bins, with pick-to-light guidance.



**3 Bulk Supplies & Kits**  
Bulk supplies, kits, and larger items are stowed neatly in configurable, locked supply drawers, with pick-to-light guidance.

### Adaptable and Scalable Configurations

NexsysADC provides a fully configurable platform to permit unique security and storage capacity for any space, daily volume, and budget. Consult with Capsa Healthcare on your precise needs.



NexsysADC is scalable for any care setting that handles controlled medications. The countertop-size NexsysADC 4T (right) holds up to 150 SKUs, and the main cabinet holds up to 400 SKUs.

### Specifications

#### Main Cabinet Model:

- Cabinet size 43"(h) 24"(d) 25"(w)
- 22" touchscreen monitor
- Capacity up to 400 SKUs (highly configurable)

#### 4T Countertop Model:

- Cabinet size 17.7"(h) 20.4"(d) 23.5"(w)
- 14" touchscreen monitor (total unit height 26")
- Capacity up to 150 SKUs (highly configurable)

### Popular configurations


Ask your Capsa representative about our common configurations to suit most facilities' needs



## Pharmacy-To-Facility Medication Management

### Cloud-Based For Continuity

NexsysADC is enabled by a robust cloud solution for superior information management and maximum uptime. From anywhere, at anytime, the NexsysADC web portal provides access to all your system management functions. All essential patient medication transactions and administrative functions are captured, such as full activity and user reporting, charge capture, and restock management.



### Pulse Interfaced With Your System

NexsysADC interfaces with your management system with our robust 2-way Pulse Interface Engine (HL7). Every step of the process is connected:

- **Patient Sync** for true patient validation reduces the chance of medication errors.
- **Patient Order Sync** only allows dispensed medications against approved orders.
- **Usage/Billing Syncs** sends usage electronically to the management system, eliminating manual billing and avoiding revenue loss.

### Uncompromising Data Security

Our NexsysADC technology meets the highest standards for cyber security and HIPAA compliance.



### Implementing NexsysADC

Capsa Healthcare makes the planning, installation and training of NexsysADC simple and straightforward. Our Support team focuses on the unique needs of clinicians.

- High-detail project management
- On-site training
- 24/7 dedicated support
- Dedicated support hotline





# AVALO<sup>®</sup> SERIES MEDICAL & PROCEDURE CARTS



## Medical & Procedure Carts

Avalo Series medical carts put everything a clinician needs at their fingertips. Avalo carts promote orderly storage, and smooth panels and seamless drawers allow simple cleaning for infection prevention. Ask about new easy-to-order convenient packages.

### Model Options



### The Avalo Medical Cart Difference

- Wide selection of accessories and organization systems
- Delivered fully assembled, no tools required
- Precision-welded steel frame with high-impact panels
- 10-year limited warranty
- See Coordinated Treatment Carts p. 23

*Available in Convenient Packages for simple ordering (left), or custom-configure.*



**Blue Emergency Cart**  
Standard Height

*Shown with Accessory Package 2*



**Red Emergency Cart**  
Intermediate Height

*Shown with Accessory Package 3*



**Procedure**



**Isolation**

### Accessories & Options



Breakaway Lock



Drawer Dividers



Sharps / Glove Holder



Oxygen Tank Holder



Defibrillator Shelf



Accessory Bridge



Bulk Storage Bins



Tilt Bin Organizers



Slide-Out Surface



Waste With Lid



Drawer Tag Holder



AC Plug Outlet Strip



Side Organizer



I.V. Pole



Cardiac Board

For additional accessories see [CapsaHealthcare.com](http://CapsaHealthcare.com)



# MOBILE COMPUTING, TELEHEALTH, VIRTUAL CARE

# Laptop & Tablet Mounting

## Point-of-Care Mobile Computing


All Capsa medication and procedure carts can be equipped for mobile computing.

We offer mounting systems for laptop, tablet and All-In-One computers, and hardware for side or back-mounting to accommodate EMR/EHR. Flexible computer positioning and adjustment save cart-top workspace.

### Laptop and Tablet Options



**AX Series Standard Laptop Tray**  
Holds up to 15" laptop




**Laptop Tray assembly, plastic, 13" x 13"**  
Holds up to 15" laptop



Lockable

**Laptop Security Tray assembly, metal**  
Holds laptops 9" - 16.5"




**Laptop Tray assembly, plastic, 13" x 17"**  
Holds up to 17" laptop



Lockable

**Laptop Plate assembly, metal, 11"x10"**  
Holds up to 15" laptop



**Tablet Holder Large**, holds 13"-16" tablet  
**Small**, holds 10"-13.5" tablet

### AX Series Cart Arms



AX Series arm & rear mount



AX Series fixed link extension



# Computer Workstations

## Documentation Carts

Capsa Healthcare offers a wide range of mobile computing workstations. Carts can be full-featured, powered or non-powered, with superior configurability to accommodate nearly any application.

### The Capsa Point-of-Care Difference

- Ideal for computerized health records and clinical documentation at the bedside or nurse station
- Fully adjustable for standing or sitting interactions
- Innovative design, sturdy construction, and aesthetic balance



### Simple Laptop Carts and Workstations

- Non-powered
- Lightweight
- Accommodates laptops, lightweight monitors, and tablets
- Simple vertical height adjustment
- Optional locking or non-locking removable drawers

*Shown: SlimCart™ (left) and Kidney Cart lightweight non-powered documentation carts*



### Powered Laptop Carts and Workstations

- Accommodates laptop, All-In-One, and tablet PCs
- Compact base for easy maneuverability
- Simple vertical height adjustment
- Large work surface
- Multiple accessories including task light, supply bins, and drawers

*Shown: Trio NP mobile computing workstation*

## Tech-Ready Medication Carts

Capsa Healthcare's tech-ready configurable cart design easily mobilizes your computing needs. Expansive capacity, durable construction, and enhanced mobility.

- Integrated slide-out keyboard and mouse drawer replaces the top 3" drawer in any standard medication cart
- Flexible computer positioning and adjustment to save cart-top workspace
- Open-platform design accommodates self-powered hardware including laptop, tablet, and All-In-One
- Multiple models accommodate a wide range of medication systems including punch card, 30-day box, and automated packaging

**Avalo® Series**  
with AX Arm  
and monitor



**M-Series**  
with AX Arm and  
All-in-One



*Also available as Vintage Encore*

# Telehealth Carts and Virtual Care Solutions

## Accommodates most tablets, iPads, monitors

Now, Capsa Healthcare delivers the widest range of point-of-care solutions that help your staff do what they do best: deliver exceptional care. Capsa's **NEW** Tryten virtual care and telehealth carts are simple and slim-profile. Tryten carts support virtual care, improved resident/patient experience, and remote interface with doctors, pharmacies, and family members.

### Virtual Care Carts

Whether you're a few floors away or across the country, virtual care carts help you bring the specialist to the resident. Tryten solutions support a wide range of medical devices and telehealth peripherals, making it easy to evaluate, assist, and supervise effectively.



### Applications

- Telehealth
- Virtual Rounding
- Patient Observation
- Virtual Nursing
- Telestroke

### Patient Experience Carts

As you shift focus to delivering improved, holistic, patient-centered care, the tools and technology you require are changing. Patient experience carts have you covered.



### Applications

- Video Remote Interpretation
- Bedside EHR
- Family Contact
- Entertainment
- Patient Registration

### Medical Device Carts

From ultrasound to telecardiology, vital sign monitoring, and beyond, we help you store, move, and use today's smaller and more powerful medical devices.



### Applications

- Vital Signs Monitoring
- Ultrasound
- Endoscope
- Cardiology



# KIRBY LESTER INVENTORY MANAGEMENT



## Inventory Management

Consider how you manage your physical inventory count, especially narcotics. Hand-counting or half bottle / full bottle estimations put your inventory integrity at risk. Capsa Healthcare provides Kirby Lester pharmacy automation to yield tighter inventory control.

### The Kirby Lester Difference

- Maintain absolute inventory accuracy (manual counting methods are 95% accurate at best)
- 14.3 seconds faster (average) than hand-counting for a 90-count
- Counting speed up to 15/second

### Multiple Applications

- Electronic C-II Log
- Perpetual inventory software
- Routine inventory counts of automated packaging robot cells
- Returns and expired medication inventory tracking
- Retail pharmacy operation, for daily Rx filling



**For The Narc Vault & Robot Cassettes: KL1 Inventory Counter**

4.6"(h) 12.5"(d) 5.6"(w); 6.9 lb



**For Retail Pharmacy Filling: KL1Plus Verify-Plus-Count Device**

14.6"(h) 12.5"(d) 8"(w); 9.6 lb

### Simple Inventory Counting: KL1



Count almost any tablet or capsule accurately & quickly



Eliminate inaccurate, inefficient manual counting

### Rx Counting/Verification: KL1Plus



For higher-level controlled medication inventory



Keep a record of Rx fills and inventory (viewable on user portal)



# Warranty, Service & Support

## Superior Products. Superior Service.

Capsa Healthcare provides a comprehensive service program focused on rapid resolution to keep your cart fleet operational, and your clinicians supported.

- Expert direct phone support available during normal business hours
- 24/7 Emergency Technical Support
- National team of Field Support Technicians
- Support via email available at [Parts.Support@CapsaHealthcare.com](mailto:Parts.Support@CapsaHealthcare.com)
- On-site technical service and support available

## Contact your Capsa Representative

for more detailed information about Capsa's warranty and service programs.



*Working in concert with your pharmacy staff and facility management team is what helps Capsa Healthcare deliver long-term customer satisfaction.*

## Contact Us

Contact us today for more information or to schedule a personal on-site demonstration.

**800.437.6633**

**International: +1.614.864.9966**

**email: [Info@CapsaHealthcare.com](mailto:Info@CapsaHealthcare.com)**

**[CapsaHealthcare.com](http://CapsaHealthcare.com)**



Certified ISO  
27001:2013



Certified ISO  
9001:2015



*Specifications, availability and components are subject to change without notice.*