

M38e



The **M38e** from Capsa Healthcare is an evolution of the most proven point-of-care computing cart in healthcare. With updated features that provide superior ergonomics, enhanced ease-of-use, and expanded storage flexibility, the M38e maximizes clinical efficiency and elevates the performance of your health IT program.

- Smaller profile glides through tight care settings
- Large work area with additional 100 square inch expandable surface feature
- Increased battery runtime for improved uptime performance
- Large CPU storage area for expanded tech organization



Ergonomic Adjustability

Simple adjustments for individual preference and ergonomic comfort



Secure

Intuitive locking and security system



Storage

Flexibility to configure workstation with storage drawers or bins

Power System	Non-Powered	AC Powered Lithium
Target Run-Time (typical configuration):	Laptop battery only	9-12 Hours/Charge
Charge Time from Empty:	Laptop charge time only	2.5 Hours
Battery Life Cycles:	N/A	2500+
Recommended Technology:	Laptop (up to 15" screen)	Laptop (folded), thin client, or All-In-One PC


Optional Accessories



Access Pack for med cups, tape and alcohol wipes and/or scanner holder



Infection Control Wipes Bracket



Printer Shelves



Forked Scanner Holder




Barcode Scanner Shelves and Scanner Holder Options



Wire Basket (standard or extra large)



Waste Bin



Zebra Printer Holder

Check with your Capsa Healthcare regional sales manager about these and more accessories specific to your chosen point-of-care solution.

Bin Configurations



Locking: FG9M3811
Non-Locking: 1874831



Locking: FG9M3833
Non-Locking: 1874836



Locking: FG9M3812
Non-Locking: 1874832



Locking: FG9M3823
Non-Locking: 1874835



Locking: FG9M3822
Non-Locking: 1874834



Locking: FG9M38DD
Non-Locking: 1874837



Locking: FG9M3813
Non-Locking: 1874833



Locking: FG9M38QQ
Non-Locking: 1874838



Pharmacy Robot Envelope Drawer

Locking Only:
1779309

Drawer Options

Flexibility to configure workstation with storage drawers in 3" or 6" depths



3" Drawer
Electronic lock or non-locking



2-3" Drawers
Electronic lock or non-locking



6" Drawer
Electronic lock or non-locking

Drawer Dimensions

3 inch: 3" w x 18.25" h
6 inch: 6" w x 18.25" h

Model Options



M38e
Document Model



M38e
XP Model



M38e
RX Model



M38e
Laptop Model

General Specifications:

- Base Size: 17" x 17"
- Weight: Configurations starting at 87.5 lbs
- Cart Height: 49.3" to 65"
- Work Surface: 17.25" W x 10.25" D
- Casters: 4" or 5" standard or directional casters
- Adjustable height, tilt/swivel monitor mount and keyboard tray options
- Enhanced lighting features to assist with navigating and working in darkened areas