

TruAir™ Scope Tracker

The TruAir™ Scope Tracker enables the efficient tracking of endoscope drying and hang time, while complying with ST91 drying guidelines. The 11" tablet with adjustable bracket easily integrates with the TruAir™ Scope Drying Cabinet.


The TruAir™ endoscope tracking software simplifies the monitoring of scope drying and hang time by identifying scope locations, storage times, and length of time before scopes expire, per your facility's infection control guidelines. Generate reports for visibility into the scope drying process and audit tracking.




Cabinet Features:

- 11-inch display tablet
- Adjustable tablet bracket for mounting on the right or left side of the cabinet
- Scope tracking software is installed on the tablet and ready for use
- Software features a color-coded system to help keep track of scope drying and hang time
- Send reports via Excel for audit tracking
- All data is saved directly on the tablet and only requires a Wi-Fi connection
- No patient or procedure information is required to use the software
- Includes one year of software upgrades and maintenance, with additional years available upon request

Scope Selection



Drying Cycle



Reporting Feature

User	Action	Status	Model
Admin	Insert	Drying	CF-HQ190L
Admin	Insert	Drying	CF-HQ190L
Chris Miller	Insert	Drying	CF-HQ190L
Chris Miller	Remove	Unallocated	CF-HQ190L
Chris Miller	Insert	Drying	CF-HQ190L
Admin	Insert	Drying	CF-HQ190L



Call 800.437.6633 for a Quote

CapsaHealthcare.com