

Tryten P1 Monitor Cart



DESIGNED FOR:


TELEHEALTH COMMUNICATION

REMOTE INTERPRETING

VIRTUAL VISITING

VIRTUAL ROUNDING

EHR WORKFLOW SUPPORT





Tryten P1 Monitor Cart

Articulating Arm

Height adjustable 7" or 14" articulating arm with cart height reaching from 44.4" to 63".

Ergonomic Handle & Cord Wrap

Ergonomic handle with slender design uniquely placed for easy push-and-position.

VESA Mount

The 75x75/100x100 VESA monitor mount supports a monitor weight of up to 9.5 lbs for the 7" arm or up to 6.5 lbs for the 14" arm, and a 110 degree tilt to position your monitor at various angles for ease of use.

Customization

Customize your cart with Branding Options and 40+ Accessories.

Best-in-Class, Durable Materials

Anodized Aluminum, Powder-Coated Aluminum, Cycology Plastic, fully wipeable and approved in operating rooms.

Concealed Cable Management

A wider column compared to the S-Series allows for cable management on both sides of the internal channel, for discreet and uninterrupted charging. Large utility box can house power cable receptacles or other power components.

Base & Casters

The 4-Wheel, larger base maximizes stability and balances the need for stability with a compact footprint. Complete with premium, quiet, non-marking, locking casters, this cart is ready to roll or stay in a locked position.

Quality & Services

5 Year Limited-Warranty, ISO 9001:2015 Certified, Tip Resistant as per IEC 60601, and Remote Expert Support.

DIMENSIONS AND SPECS

- Column Articulating Arm: 44.4" to 63" H
- Large Utility Box, Inside fit, 8.6" x 3.9" x 2.0"
- Large Utility Box, Outside, 9.2" x 5.5" x 2.8"
- Base footprint: 18.7" to 20.6"
- Monitor weight limit:
 - (7" arm): 9.5lbs
 - (14" arm): 6.5lbs

WHAT'S INCLUDED

- P-Series Poll with Articulating Arm and VESA Bracket
- Base with Locking Casters
- Large Utility Box
- Handle and Cord Wrap
- Hardware Pack
- Assembly Instructions